

Score Medical Line stools

The Score Medical Line stools are available in two shapes: a traditional complete round seat and an ergo-shape seat with recessed edges to prevent pinching of the thighs.



Score Medical 6211 Round seat



Score Medical 6311 Ergo-shape seat



Technical details:

- Height adjustment seat standard gas lift MID31 : 56 – 75 cm
- Height adjustment seat high gas lift MID38 : 63 – 88 cm
- Height adjustment seat low gas lift MID25 : 47 – 61 cm
- Dimensions seat (diameter) : Ø 40 cm
- Thickness seat : 7 cm
- Dimensions lumbar support (height x width) : 25 cm x 39 cm
- Mechanism : Tilt
- Inclination adjustment seat : +6° / -6°
- Inclination adjustment lumbar support : +11° / -21°
- Diameter base : 54 cm / optional 50 cm
- Diameter castors : 65 mm
- Maximum weight use : 120 kg
- Guarantee : 5 years

Score Medical 6260 Round seat



Score Medical 6360 Ergo-shape seat



Technical details:

- Height adjustment seat standard gas lift MID31 : 54 – 74 cm
- Height adjustment seat high gas lift MID38 : 61 – 86 cm
- Height adjustment seat low gas lift MID25 : 46 – 61 cm
- Dimensions seat (diameter) : Ø 40 cm
- Thickness seat : 7 cm
- Mechanisms : Fixed
- Diameter base : 67 cm / optional 60 cm
- Diameter castors : 65 mm
- Maximum weight use : 120 kg
- Guarantee : 5 years

Options:

- Height adjustment by foot control
- Foot ring Ø 45 cm

Upholstery: Stamskin or PUxx

